

CPAC Autex Medical X-Ray Chemistry

Autex Developer starter

For use with the X-Ray developer to make a working tank solution in any processor. With 1 Liter concentrate you can make 40 Liter of working tank solution (=25 ml/liter).

Catalog Number _____ **Mix Size** _____ **Shipping Weight**

20DE2010 _____ 6 x 1Lit Concentrate _____ 7,50 kg

Pallet info: 150 kits per pallet

CPAC Autex LC Developer

This developer sets the quality standard for general radiography. Autex® LC Developer and Replenisher is a three-part, liquid concentrate. It processes all types and brands of medical X-ray film in automatic roller transport processors with processing cycles of 90 seconds, to over 3 minutes 30 seconds. Autex LC is available in convenient 2 x 20 BiPak® (40)liter mixes

Features:

- Suitable for all brands of medical X-ray film
- Clean working
- Cost effective
- Balanced to provide low fog, high speed and good toe contrast
- Excellent latitude
- Formulated for long shelf life
- High D-max

Catalog Number _____ **Mix Size** _____ **Shipping Weight**

20DE2200 _____ 2 x 20 L _____ 14,20 kg

Pallet info: 68 Kits per pallet

Concentrate specifications:

Part A is a 5 Lit concentrate

Part B is a 500ml concentrate

Add 14,5 Lit water to this to make 20Liters of Ready to use developer

Each Kit contains 2 Part A & B

CPAC Autex LC Fixer

Autex® automatic LC X-ray fixer and Replenisher is two-part, all-liquid hardening fixer designed to provide superior fixing and hardening for all medical and industrial X-ray films. It is especially formulated for high temperature automatic processors, but will perform equally well as a hand tank fixer. Its low-odor formulation makes it great for smaller or older darkroom environments. This low-odor formula does not use a scent to mask the odor. Our unique formula actually reduces the fumes that cause fixer odor. Autex® Automatic LC X-ray Fixer and Replenisher is available in convenient 2 x 20 BiPak® liter (40 liter) mixes

Features:

- Suitable for types and brands of medical and industrial X-ray film
- Clean working
- Non-corrosive
- Designed for high temperature automatic processing
- Non-sludging

Catalog Number _____ **Mix Size** _____ **Shipping Weight**

20FX2400 _____ 2 x 20 L _____ 12.00 kg

Pallet info: 72 Kits per pallet

Concentrate specifications:

Part A is a 4 Liter concentrate

Part B is a 500ml concentrate

Add 15,5 Lit water to this to make 20Liters of Ready to use fixer

Each Kit contains 2 Part A & B.