

Medical X-Ray product range

Autex Developer starter

For use with the X-Ray developer to make a working tank solution in any processor. With 1 Liter concentrate you can make 40 Liter of working tank solution (=25 ml/liter).

Catalog Number	Mix Size	Shipping Weight
20DE2010	6 x 1Lit Concentrate	7,50 kg

CPAC Autex Manual Developer

This developer sets the quality standard for general radiography. Autex® Manual Developer and Replenisher is a one-part, liquid concentrate. It processes all types and brands of X-ray films. It is suitable for processing cycles of 90 seconds to over 6 minutes.

It is designed to ensure maximum life/stability, tank cleanliness and optimal sensitometric film performance. The Autex Manual Developer is formulated for superior manual processing. This premium developer is suitable for use with all types and brands of medical, dental, industrial and cinefluorograph films.

Autex Manual Developer is available in convenient 6 x 5L Packs (30) liter mixes.

Features:

- Suitable for all brands of X-ray film
- Low fog with balanced speed and contrast for maximum detail
- Clean working
- Cost effective
- Consistent image quality
- Minimum replenishment
- High capacity
- Excellent latitude
- Formulated for long shelf life

Catalog Number	Mix Size	Shipping Weight
20DE2001	6 x 5 L	10.50 kg